

THE INDUSTRY'S MOST EFFICIENT LARGE-CAPACITY TABLE-TOP STERILIZER

RH-Pro11

HIGH-VELOCITY HOT AIR (HVHA) STERILIZER

NO STEAM
NO DRYING
NO CORROSION
NO MAINTENANCE



MADE IN THE
USA

The RH-Pro11 Advantage

- Waterless Sterilization with No Chemicals
- Complete 12-Minute Wrapped Cycle
- Uses 85% Less Energy than Steam
- Accommodates Small and Large Cassettes
- Up to Six (6) Large Instrument Trays Per Cycle
- Documents Cycles for Easy Log Retrieval
- Self-Locking Door Insures Cycle Completion

Specifications

Unit Size (OD)

22.5" W x 22.5" D x 19.5" H
572mm x 572mm x 495mm

Chamber Size (ID)

11"W x 17.75"D x 11.75"H
279mm x 433mm x 299mm

Chamber Capacity

2294 cubic inches (10 gal/38 liters)

Instrument Tray Size (ID)

9"W x 15"D x 1"H
229mm x 381mm x 28mm

Unit Weight: 75 lbs (34kg)

Electrical

110-120V, 60Hz, 1400 Watts warm up,
300 Watts operating
220-240V, 50Hz, 1400 Watts warm up,
300 Watts operating

Sterilization Temperature

375° Fahrenheit (190° Celsius)

Warranty

Three-Year (parts and labor)

Approval/Clearance

FDA 510(k) K881371
TUV Electrical/EMI Safety,
CE Directive Compliance

Sterilization Times

Cycle 1: Six (6) minutes for unwrapped instruments

Cycle 2: Eight (8) minutes for hand pieces and medical drills

Cycle 3: 12 minutes for pouched, wrapped instruments, trays and cassettes

Included with Sterilizer:

Instructions & Care Card
User Manual
Four (4) Instrument Trays
One (1) Pair Tray Removal Safety Gloves

Equipment Options:

Cassette & Pouch Racks

CPAC
ENVIRONMENTAL SOLUTIONS

RH-Pro11
HIGH-VELOCITY HOT AIR STERILIZER

Store.Sterilizers.com

(866) 309-8654